



CENTER FOR BIOMEDICAL IMAGING IN ONCOLOGY

Lurie Family Imaging Center



Dana-Farber
Cancer Institute



HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL

DATA PREPARATION: Instruction for machine learning data preparation

STEP 1

- Run ClinicalVolumes tumor segmentation. Save segmentation of the tumor into a “seg” folder. one “seg” folder for each animal

STEP 2

- **RGB image creation (Python)**
 - 1) Open “Dicom_to_tiff/1_RGBcreator/Python” folder
 - 2) Inside the “data_dir” folder, create a new directory (Note: create one directory for each animal)
Name each directory similarly to the following: LFIC_MR_0001_20200101
The name must have the same number of letters and underscore as the suggested. The name will be used to name the RGB images
 - 3) Move or copy the dicom files into the corresponding folder inside the “data_dir” directory
 - 4) Open “RGB_creator_multidir” in a Python editor
Select image options as requested
 - 5) Run “RGB_creator.py” from the terminal with the command:

```
python RGB_creator_multidir.py
```
- The program will convert all the dicom files inserted in the different folders in .tif 256x256x3

STEP 3

- **8bits mask image creation (MATLAB)**
 - 1) Copy a single “seg” folder (folder with the segmentation dicom files) into “Dicom_to_tiff/1_RGBcreator/Matlab” folder
 - 2) Open Gray_scale_and_8uint_image.m file in MATLAB and follow the instruction inside the file.
(no need to clean the folders from old data)
 - 3) Segmented masks will be created and found in the “Images_GrayS_8bits” folder

- The program must be run for one animal (one dicom dataset) at the time

Note: If needed to **rename**, **invert the order** or **rotate** the .tif image files:

1. open the "Dicom_to_tiff/ 3_IMAGE_editing" folder
2. move the tif images into the corresponding folder (The name of the folder will be used to name or rename the images)
3. Open "IMAGE_editing_multidir.py" in a Python editor
Select image options as requested
4. Run "IMAGE_editing_multidir.py" from the terminal with the command: `python IMAGE_editing_multidir.py`

STEP 4

- Move the tif image files and mask files in a folder XX. One folder per animal



CENTER FOR BIOMEDICAL IMAGING IN ONCOLOGY

Lurie Family Imaging Center



Dana-Farber
Cancer Institute



HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL

Instruction for machine learning data preparation (Skull_stripping)

STEP 1

- Run ClinicalVolumes tumor segmentation. Save segmentation of the tumor into a “seg” folder.
- Run ClinicalVolumes full brain segmentation. Save segmentation of the brain into a “seg_skull_stripping” folder.

STEP 2

- **RGB image creation _ Skull_stripping_mask creation (MATLAB)**
 - 1) Move the “seg_skull_stripping” folder into
“Dicom_to_tiff/2_RGBcreator_Skull_stripping/Matlab” folder
 - 2) Open Skull_stripping_mat file in MATLAB and follow the instruction reported in the file. (no need to clean the folders from old data)
 - 3) Segmented brains .tif files will be found in the “Skull_stripping_mask” folder
- The program must be run for one animal (one dicom dataset) at the time
- **RGB image creation _ Skull_stripping_mask creation (Python)**
 - 4) Open “Dicom_to_tiff/2_RGBcreator_Skull_stripping /Python” folder
 - 5) Inside the “data_dir” folder create a new directory (Note: create one directory for each animal)
Name each directory similarly to the following: LFIC_MR_0001_20200101
The name must have the same number of letters and underscore as the suggested. The name will be used to name the RGB images
 - 6) Inside the “Skull_stripping_mask” folder create a directory of the same name as in (5)
 - 7) Move or copy the dicom files into the corresponding “data_dir” folder
 - 8) Move or copy the .tif skull stripping masks files (created in STEP 2) into the corresponding
“Skull_stripping_mask” folder

- 9) Open "RGB_creator_multidir_stripping.py" in a Python editor
Select image options as requested
- 10) Run "RGB_creator_multidir_stripping.py" from the terminal with the command:
python RGB_creator_multidir_stripping.py
- The program will convert all the dicom files inserted in the different folders in .tif 256x256x3

STEP 3

- **8bits mask image creation (MATLAB)**
 - 4) Copy the "seg" folder into "Dicom_to_tiff/1_RGBcreator/Matlab" folder
 - 5) Open Gray_scale_and_8uint_image.m file in MATLAB and follow the instruction inside the file.
(no need to clean the folders from old data)
 - 6) Images will be created in the "Images_GrayS_8bits" folder
- The program must be run for one animal (one dicom dataset) at the time

Note: If needed to **rename**, **invert the order** or **rotate** the tif image files:

5. open the "Dicom_to_tiff/ 3_IMAGE_editing" folder
6. move the tif images into the corresponding folder (The name of the folder will be used to name or rename the images)
7. Open "IMAGE_editing_multidir.py" in a Python editor
Insert the directory paths as requested
Select image options as requested
8. Run "IMAGE_editing_multidir.py" from the terminal with the command: *python IMAGE_editing_multidir.py*

STEP 4

- Move the .tif image files and .tif mask files in a folder XX. One folder per animal